

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22604-1K	2026-08-31
Low	22604-11	2026-08-31
Normal	22604-12	2026-08-31
High	22604-13	2026-08-31

Open vial stability 14 days <sup>1</sup>

## SWELAB

CONTROL		L	
LOT		22604-11	
Swelab Alfa/Alfa Plus		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.23	0.12
MCV	fL	67.1	5.0
RDW	%	14.0	4.0
Hct	%	15.0	2.5
Plt	10 <sup>9</sup> /L	81	15
MPV	fL	11.6	1.5
WBC	10 <sup>9</sup> /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 <sup>9</sup> /L	1.6	0.4
Mid#	10 <sup>9</sup> /L	0.1	0.3
Gran#	10 <sup>9</sup> /L	1.7	0.8
Lymph	%	46.7	7.0
Mid	%	3.5	5.0
Gran	%	49.8	7.0

CONTROL		N	
LOT		22604-12	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		4.20	0.18
		85.1	5.0
		12.1	4.0
		35.7	3.0
		238	30
		11.8	1.5
		8.4	0.6
		12.3	0.4
		29.3	1.5
		34.5	3.0
		4.2	0.6
		0.4	0.5
		3.8	1.0
		50.3	5.0
		4.2	5.0
		45.5	5.0

CONTROL		H	
LOT		22604-13	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		5.00	0.22
		93.9	5.0
		12.0	4.0
		47.0	3.5
		552	60
		12.3	1.5
		20.6	1.8
		15.9	0.5
		31.8	1.5
		33.8	3.0
		13.2	1.5
		0.8	1.0
		6.6	2.5
		64.5	5.0
		3.7	5.0
		31.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonu stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
  - Sredstvo
  - Průměr
  - Middelværdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Średnia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Rasjon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Orseg
  - Rozsah
  - Območje
  - Intervalo
  - Intervall
  - Aralığı



[www.boule.com](http://www.boule.com)

 **Boule**



Boule Medical AB  
Fagerstagatan 7  
SE-163 53 Spånga,  
Sweden

25121-8

## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa  
 Lot: 22604-11  
 Prod.: Con-Diff Low  
 ID: 2260411+  
 Exp.: 31-Aug-2026

Instr.: Swelab Alfa  
 Lot: 22604-12  
 Prod.: Con-Diff Norm  
 ID: 2260412+  
 Exp.: 31-Aug-2026

Instr.: Swelab Alfa  
 Lot: 22604-13  
 Prod.: Con-Diff High  
 ID: 2260413+  
 Exp.: 31-Aug-2026

